

accelerate productivity with unequalled durability



Thermo Scientific[™] Fiberlite[™] rotors maximize centrifuge performance with versatility, speed and a robust corrosion-free design

Improved ergonomics and productivity

Lightweight design

Large metal centrifuge rotors often present a unique lifting hazard in the laboratory due to their weight and awkward shape. Lightweight Fiberlite rotors – up to 60% less weight than metallic rotors¹ – feature improved ergonomics for a safer work environment and minimize the risk of damage to centrifuge equipment.

Additionally, these lightweight properties result in faster acceleration/deceleration rates for shorter run times.

Unequalled durability and cleaning convenience

Corrosion and fatigue resistance

Traditionally, the primary cause of rotor failure is from damage to metal surfaces due to moisture, chemicals or alkaline solutions that weaken the metal rotor's structural integrity. Carbon fiber composite rotors are corrosion-resistant, eliminating this ever-present hazard, and are safe to use with most mild laboratory detergents and solutions, ensuring easy rotor care and maintenance.

Substantial load or stress, as a result of high rotational speeds and repeat cycles, can also threaten metal rotor structure by causing it to stretch and change in size, limiting rotor life or leading to failure. Fiberlite rotors are fatigue-resistant, eliminating this threat.

Exceptional value within your reach

15-year warranty² in all centrifuges

Unlike the limited lifetime of metal rotors due to potential failure risks, Fiberlite carbon fiber rotors are backed by the most comprehensive warranty² coverage available.

Unique repairability

In contrast to traditional metal rotors, Fiberlite carbon fiber rotors are repairable if damaged.

Superior insulation

Carbon fiber material possesses naturally insulating properties, which helps to maintain sample temperature integrity.

Weight comparison of fully loaded 6 x 250 mL capacity floor model rotors¹



ROTOR WEIGHT KG/LBS

¹ Based on a comparison with manufacturers' published specifications.

² Subject to Thermo Fisher Scientific's standard limited warranty. See thermoscientific.com or your sales representative for details.

Figure 1: Weight savings with carbon fiber rotors.

Fiberlite rotors are available for a wide range of processing needs -

ADME/Toxicology **Bioproduction Blood Banking Cell Biology Cell Culture Cell and Tissue Analysis Forensic Analysis** Formulation Immunology

Microbial Testing

Microbiology

Nucleic Acid Research

Nutritional and Dietary Concerns

Pathology

Pharmaceutical QC and Production

Protein Analysis, Isolation and Expression

RNAi and Gene Regulation

Stem Cell

Water and Waste Water Analysis and Water Pollution Analysis

Acceleration and deceleration rate comparison of 6 x 250 mL capacity floor model rotors¹

specific warranty coverage for each rotor.

manufacturers' published specifications.

³ Warranty coverage may vary by rotor. Please refer to manufacturer for

⁴ Average warranty periods were calculated based on industry average of years an aluminum or titanium rotor may be covered under warranty per



				with Fiberlite carbon fiber rotors*											
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	
Fib	erlite	Carl	on	Fiber	r Rot	tor									
Alu	minuı	n Ro	tor ³												
Tita	nium	Roto	r3												

Average warranty periods for metal rotors compared

WARRANTIED ROTOR LIFE

OPERATION WITH PERIODIC INSPECTIONS

Figure 2: Time savings with carbon fiber rotors.

Figure 3: Warranty with carbon fiber rotors.

Best-in-class Thermo Scientific Fiberlite



rotor portfolio

Seamless integration

From benchtop instruments to advanced floor models, Thermo Scientific centrifuge systems deliver outstanding performance and reliability in the lab. We provide an integrated solution of rotors, equipment, and accessories, offering exceptional value and best-in-class features including:

- innovation and technical design
- high throughput and speed
- operator, sample and system safety
- · operational longevity of your system

Superior sample containment

- In the event of a tube or bottle failure, a volume of fluid is contained inside the rotor in a liquid containment annulus, preventing biohazardous samples from escaping; available on select rotors.
- To enhance containment of biohazardous samples, rotors certified by the Public Health Laboratory Service, Microbiology Services, Porton Down, UK are noted by
- Lids for rotors featuring Thermo Scientific[™] Auto-Lock[™] rotor exchange enable rotors to remain sealed while being carried to a biocontainment hood for sample retrieval; available on select rotors.



Thermo Scientific Fiberlite LEX Rotor Series

Introducing the latest innovation in Thermo Scientific Fiberlite carbon fiber rotor technology.

The next generation of high capacity Fiberlite rotors, the Fiberlite LEX rotor series, further advances the current carbon fiber design, combining even lower mass with low kinetic energy to deliver superior ergonomics with outstanding performance and safety.

Fiberlite	F9-6x1000 LEX
	F10-4x1000 LEX
	F12-6x500 LEX
	F20-12x50 LEX

Superior ergonomics

Fiberlite LEX rotors take the lightweight design of carbon fiber to a whole new level; these rotors are the lightest of their kind¹, with improved ergonomics for everyday ease of handling.

Exceptional performance

The Fiberlite LEX rotor series provides outstanding RCF performance for enhanced productivity – up to 24,471 x g with the 6 x 500 mL (3-Liter volume) Fiberlite LEX rotor and up to 17,568 x g with the 6 x 1000 mL (6-Liter volume) Fiberlite LEX rotor.

Weight comparison of 6 x 500 mL and 6 x 1000 mL capacity floor model rotors¹



Figure 4: Lower weight advantage of Fiberlite LEX rotors.



Enhanced safety

From sample protection with advanced sealing properties, to safety of equipment and lab personnel with the rotor's lifting handle, Fiberlite LEX rotors are the top choice for a safe work environment.

In today's biomedical and microbiological laboratories, containment of biological agents and infectious substances are an essential element in maintaining a safe environment. Fiberlite LEX rotors provide multiple levels of protection to enhance biosafety without compromising functionality or convenience.

- Biocontainment Tested: Rotors certified by the Public 1 Health Laboratory Service, Microbiology Services, Porton Down, UK are noted by
- 2 Liquid Containment Annulus: In the event of a bottle failure, a volume of fluid is contained inside the rotor, preventing biohazardous samples from escaping.
- 3 Auto-Lock Rotor Exchange with Thermo Scientific[™] Auto-ID[™] Rotor Identification: Simplifies run set-up and eliminates the worry of overspeeding or rotor accidents.

Lower kinetic energy resulting from the lightweight design, enhances equipment performance and safety of work environment.



¹ Based on a comparison with manufacturers' published specifications

Superspeed Rotors

With volumes ranging from 1.5 mL to 6 Liters, a full range of Fiberlite carbon fiber rotors is available for superspeed floor model centrifuges, facilitating applications spanning pharmaceutical, biotechnology and academic research.



the built-in lifting handle and liquid containment annulus (available on select rotors).

Fiberlite rotor model nomenclature



High capacity and seamless compatibility

Fiberlite F9-6x1000 LEX F10-4x1000 LEX F12-6x500 LEX F14-6x250y

- Simplify preparation by loading tubes directly into Fiberlite rotors, eliminating multi-piece canister assemblies, which can be misplaced or damaged.
- Work seamlessly with Thermo Scientific[™] Nalgene[™] bottles, including the high performance 1-Liter widemouth polypropylene and polycarbonate centrifuge bottles that process a full Liter at maximum speeds (20,584 x g) with leakproof assembly.

Enhanced ergonomics

- Lightweight design allows easy rotor transport in and out of the centrifuge.
- Installation or exchange of rotors requires less force

 especially with lifting handle on select models –
 reducing risk of injury.

Conical tube efficiency

Fiberlite | F14-14x50cy

- Spin 14 x 50 mL conical tubes at maximum rotor speed (33,700 x g) without tube damage.
- Process 15 mL conicals with available adapters for flexibility.

Small-volume protocol support

Fiberlite F20-12x50 LEX F21-8x50y F23-48x1.5

• Small-volume pelleting and microtubes ranging from 1.5 to 50 mL at RCFs up to 57,300 x g.

⁶ Actual maximum rotor speed may vary depending on centrifuge.

7 Actual fill volumes may vary from nominal volume.

Fiberlite Rotors for the **NEW! Thermo Scientific[™] Sorvall[™] LYNX** Superspeed Centrifuge Series

Rotor innovations shorten run set-up time while providing peace-of-mind that the rotor is secure.



Figure 7: Auto-Lock rotor exchange: Secure, trouble-free rotor installation and removal in as little as 3 seconds.





Figure 9: Speed handle on rotor lids: Makes tightening the lid safer while also simplifying lid removal.

Figure 8: Auto-ID instant rotor identification: Improves safety, saves times, and protects the integrity of your samples.

Auto-Lock rotor exchange

Secure, push-button rotor exchange in as little as 3 seconds delivers:

- Improved safety and confidence that the rotor is automatically and securely locked and will not loosen during a run
- Trouble-free rotor installation and removal
 - 1 No tools are required
 - 2 The rotor locks itself to the centrifuge, eliminating the need for hand-tightening
- Flexibility to quickly change rotors and applications, matching the needs of your laboratory – today and in the future

Auto-ID instant rotor identification

Immediate identification of a rotor when secured in the centrifuge chamber, with rotor specifications automatically loaded into the centrifuge parameters.

- Shortens run set-up time by eliminating the need to find and set rotor codes
- Eliminates over-speed risk, reduces error messages, and improves centrifuge, sample and operator safety

Speed handle on rotor lids

- Accelerates and simplifies rotor lid tightening, ensuring lid is properly attached
- Easier and safer lifting and carrying of rotors, further enhanced with the lightweight design

innovative rotor convenience

Conical Tubes

Complete workflow in disposable conical tubes

Fiberlite F13-14x50cy F14-14x50cy F15-8x50cy

- Run samples in inexpensive disposable conical tubes, protecting from contamination and reducing sample transfers and non-productive tasks, such as autoclaving.
- Reduce processing times by spinning at maximum speeds up to 33,700 x g⁸ without risk of tube damage.
- Clarify crude lysates for plasmid DNA preps from Qiagen[™] Maxi and Midi Prep protocols.

Figure 10: Through exclusive technology, Fiberlite rotor cavities are molded to the exact shape of many disposable conical tubes for maximum support; 50 mL conical tube shown here. In addition, a cap support is designed to relieve high g-forces.

Spin sample in one tube until it's ready to store.



Figure 11: Support preparative centrifugation in a single conical tube for time and cost efficiencies and waste reduction.

Ultraspeed Rotors

From proteomics and cell clarification to nucleic acid preparation, the superior design and manufacturing of Fiberlite ultraspeed rotors deliver high performance, eliminating corrosion and the need for derating or reducing speed over the rotor lifespan.

Large volume processing

Fiberlite F37L-8x100

- Realize 33% more capacity¹ with two additional tube cavities for high volume separations.
- Achieve forces of up to 182,460 x g for time savings on separations of subcellular organelles or concentration of viruses.
- Collect or purify small macro molecular species including enzymes, antibodies and proteins from standard culture flasks up to 500 mL in a single run.

Remarkable sample throughput of microtubes

Fiberlite F50L-24x1.5

- Provide full tube support at RCF of 280,000 x g for sharp and efficient pelleting of microparticles in high performance microtubes.
- Run partial filled tubes, as low as 0.2 mL, at maximum speed for extended times without excessive tube crazing or sample loss.
- Experience multifunctional use for preparative analysis with ultracentrifuge systems.

¹ Based on a comparison with manufacturers' published specifications.

Figure 12: Fiberlite ultraspeed rotors (counterclockwise from top right): F37L-8x100 (37,000 rpm; 182,460 x g); F50L-8x39 (50,000 rpm; 266,280 x g); F50L-24x1.5 (50,000 rpm; 280,000 x g); F65L-6x13.5 (65,000 rpm; 324,140 x g).



Large Capacity Rotors

Fiberlite large capacity rotors are ideal for batch bioprocessing of bacteria or yeast and clinical samples.

Figure 13: Fiberlite F8-6x1000y rotor (8,500 rpm; 15,900 x g).

Thermo F8-5x1000y

Benchtop Rotors

Choose a Fiberlite benchtop rotor solution for high speed applications including PCR post-reaction cleanup, cell culture, plasma and general purpose separations, DNA sample preparation, subcellular fractionation and protein identification.

Accelerated applications

Fiberlite | F14-6x250LE | F15-6x100y

 Achieve outstanding g-force without compromising capacity – 250 mL up to 18,500 x g; 100 mL up to 24,500 x g – allowing more applications to be done on the benchtop.

Conical tube efficiency

Fiberlite | F13-14x50cy | F15-8x50cy

- Provide generous 14- or 8-place 50 mL conical tube capacity, and g-forces up to 24,446 x g for sample preparation without the need to change tubes.
- Process 15 mL conical tubes with optional adapters for flexibility.

Micro-volume protocol support

Fiberlite | F21-48x1.5

- Run up to 48 tubes at over 25,000 x g, doubling the capacity of conventional rotors and reducing processing by half.
- Provide ultimate user convenience with non-corroding, dual-row configuration.

Outstanding microplate processing

Fiberlite | H3-LV

- Experience exceptional capacity of 28 standard plates or 8 deep-well plates per run with g-forces up to 2,740 x g.
- Ideal for pelleting cells and cellular debris, protein precipitation and collecting physiological fluids for diagnostic testing.





Figure 14: Easy and secure push-button Auto-Lock rotor exchange in as little as 3 seconds for application versatility and cleaning convenience.

Specifications/Ordering Information

Rotors		C Sample		Cat. No.	Related Centrifuge	Max Speed	Max RCF
Rotors		Containment		Gal. NO.	Thermo Scientific	(rpm)	(x g)
Sorvall LYNX	Superspe	ed	Rotors w	ith Auto-L	ock		
	Fiberlite F9-6x1000 LEX		\ ∲	096-061075	Sorvall LYNX 6000	9,000	17,568
	Fiberlite F10-4x1000 LEX		\ ∲	096-041075	Sorvall LYNX 6000, 4000	10,500	20,584
	Fiberlite F12-6x500 LEX		\ ∲	096-062375	Sorvall LYNX 6000, 4000	12,000	24,471
	Fiberlite F14-6x250y		\ €	096-062075	Sorvall LYNX 6000, 4000	14,000	30,240
	Fiberlite				Sorvall LYNX 6000	14,000	33,746
	F14-14x50cy	•	⊗	096-145075	Sorvall LYNX 4000	13,000	29,097
	Fiberlite				Sorvall LYNX 6000	20,000	51,428
	Fiberlite F20-12x50 LEX		⊗	096-124375	Sorvall LYNX 4000	18,000	41,657
	Fiberlite				Sorvall LYNX 6000	20,000	47,850
	F21-8x50y		⊗	096-084275	Sorvall LYNX 4000	18,000	38,759
	Fiberlite				Sorvall LYNX 6000	23,000	57,368
	Fiberlite F23-48x1.5			096-484075	Sorvall LYNX 4000	18,500	37,116
C = Conical Tubes Biocontainment certification by the Public Health Laboratory Service,							

Ÿ Microbiology Services, Porton Down, UK.

Specifications/Ordering Information

Potore		c	Sample	Cat. No.	Related Centrif	Max Speed	Max	
Rotors		C	Containment		Thermo Scientific	Beckman™	(rpm)	RCF (x g)
Superspeed	Rotors							
8	Fiberlite F8-6x1000y			76641	SorvalI [™] Evolution [™] RC Series		8,500	15,800
					Sorvall [™] RC 6 [™] Plus		9,500	16,880
12.1	Fiberlite		- 🍫	096-041053	Sorvall Evolution RC Series		9,000	15,150
	F10-4x1000 LEX				Sorvall RC-5, RC-2 Series		7,000	9,160
A	Fiberlite F12-6x500 LEX		÷	096-062185	Sorvall RC 6 Plus, Evolution RC Series		12,000	24,500
					Sorvall RC-5, RC-2 Series		10,000	17,000
	Fiberlite F10-6x500y			096-062114		J2, Avanti [™] Series ⁹	10,000	17,700
	Fiberlite F14-6x250y		÷	78500	Sorvall RC 6 Plus, Evolution RC, RC-6, RC-5, RC-2 Series		14,000	30,100
(EEE)	Fileerlite			46922	Sorvall RC 6 Plus, RC-5, RC-2 Series		13,000	29,000
	Fiberlite F13-14x50cy	•		096-145011		J2, Avanti Series ⁹	14,000	33,600
	Fiberlite F20-6x100			096-064025	Sorvall RC 6 Plus, RC-5, RC-2 Series		20,000	43,900
	Fiberlite F21-8x50y		÷	46923	Sorvall RC 6 Plus, RC-5, RC-2 Series		20,000	47,500
	Fiberlite F21-48x1.5			096-484020	Sorvall RC 6 Plus, RC-5, RC-2 Series		20,000	43,500
		Capical Tubas	Discussion of the difference of the discussion o	Dublic Health Laboratory				

⁹ Except the Avanti J-HC.

 $\mathbf{C} = \text{Conical Tubes}$

Biocontainment certification by the Public Health Laboratory Service, Microbiology Services, Porton Down, UK.

Rotors			Sample	Cat. No.	Related Centrifu	Max Speed	Max RCF		
		C	Containment	Gal. NO.	Thermo Scientific	Beckman	(rpm)	(x g)	
Ultraspeed Re	otors								
8	Fiberlite F37L-8x100			096-08056	Sorvall WX Series	L8 Series ¹⁰	37,000	182,460	
	Fiberlite F50L-8x39			096-087051	Sorvall WX Series	L8 Series ¹⁰	50,000	266,280	
	Fiberlite F65L-6x13.5			096-067135	Sorvall WX Series	L8 Series ¹⁰	65,000	324,140	
	Fiberlite F50L-24x1.5			096-247028	Sorvall WX Series	L8 Series ¹⁰	50,000	280,000	
Large Capacity Rotors									
8	Fiberlite F8-6x1000y			096-061137	Sorvall RC BIOS		8,500	15,900	
	Fiberlite H3-LV			096-028016	Sorvall RC 3B, RC 3C Series		3,200	1,940	
Benchtop Rot	tors					·			
Fiberlite F14-6x250 LE			*	75003662	Sorvall [™] Legend [™] XT, Heraeus [™] Multifuge [™] X3, SL 40F Series		10,000/ 11,000 ¹¹	15,317/ 18,533 ¹¹	
	Fiberlite				Sorvall Legend X1, Sorvall Legend X X1, Heraeus Multifuge X3, SL 40F S	T, Heraeus Multifuge	15,000	24,652	
	F15-6x100y		*	75003698	Sorvall ST 16, Sorvall ST 40, Heraeu SL 16, SL 40 Series	13,000	18,516		
					Sorvall Legend X1, Heraeus Multifuge	X1 Series	8,500	12,359	
	Fiberlite	-		75003661	Sorvall Legend XT, Heraeus Multifuge	X3, SL 40F Series	9,250/ 10,000 ¹¹	14,636/ 17,105 ¹¹	
	F13-14x50cy			75006526	Sorvall Legend T, Heraeus Multifuge	9,250/ 10,000 ¹¹	14,636/ 17,105 ¹¹		
				75003663	Sorvall Legend X1, Sorvall Legend XT, Heraeus Multifuge X1, Heraeus Multifuge X3 Series		14,500	24,446	
	Fiberlite F15-8x50cy	-		75006516	Sorvall Legend T, Heraeus Multifuge	12,000/ 14,500 ¹¹	16,741 / 24,446 ¹¹		
Fiberlite F21-48x1.5/2.0			�	75003664	Sorvall Legend X1, Sorvall Legend X3, Sorvall ST 40, Heraeus Multifuge X1, Heraeus Multifuge X3, Heraeus Megafuge 40, SL 40 Series		15,200	25,055	
				75006527	Sorvall Legend T, Heraeus Multifuge	3 Series	15,000	24,400	
	Fiberlite H3-LV			75003665	Sorvall Legend XT, Heraeus Multifug SL 40F Series	je X3,	3,600	2,738	

9 Recommended for L8M and prior models.
 10 With 230 V centrifuge

 $\bm{C}=\text{Conical Tubes}$

Biocontainment certification by the Public Health Laboratory Service, Microbiology Services, Porton Down, UK.

Perfect Fit

Thermo Scientific Fiberlite rotors with Nalgene bottles and tubes **bring together quality and performance.**

Select Fiberlite rotors come complete with an initial set of Nalgene products.

Nalgene Bottles and Tubes	Nominal Capacity ⁷ per Cavity	Description	Cat. No.	Fiberlite Rotor
	1L	Nalgene Wide-mouth Superspeed Bottle, PC; SCA, PP	3140-1002	F9-6x1000 LEX
		Nalgene Wide-mouth Superspeed Bottle, PPCO; SCA, PP	3141-1002	F10-4x1000 LEX
	500 mL	Nalgene Wide-mouth Superspeed Bottle, PC; SCA, PP	3140-0500	F12-6x500 LEX
		Nalgene Wide-mouth Superspeed Bottle, PPCO; SCA, PP	3141-0500	
	250 mL	Nalgene Bottle, PC; SCA, PP	3140-0250	F14-6x250y
	250 IIIL	Nalgene Bottle, PPCO; SCA, PP	3141-0250	F14-6x250 LE
	50 ml	Nalgene Oak Ridge Tube, PC; SCA, PP	3138-0050	F21-8x50y
U	50 mL	Nalgene Oak Ridge Tube, PPCO; SCA, PP	3139-0050	F20-12x50 LEX

PC = Polycarbonate PPC0 = Polypropylene copolymer SCA = Screw closure assembly PP = Polypropylene

7 Actual fill volumes may vary from nominal volume.

Optimize the performance of your centrifuge

It's simple. From 1 L bottles, to 15 and 50 mL conical tubes, to microplates and tissue culture flasks, the versatile selection of **Thermo Scientific Nalgene and Nunc centrifugation products** work seamlessly with your complete centrifuge and rotor system, bringing together quality and performance.



Thermo Scientific Fiberlite Rotor Adapters and Accessories

	Rotor Volume ⁷ Description	No. of Vessels per Adapter	Cat. No.	Rotor Volume ⁷ Description	No. of Vessels per Adapter	Cat. No.
	•					
39 ml	. Ultraspeed			500 mL		
1:	3.5 mL Tube	1	010-1142	250 mL Conical Tube	1	010-1135
				250 mL Oak Ridge Tube	1	010-0151
50 ml				175 mL Nalgene Conical Bo	ottle 1	010-0152
_	0 mL Oak Ridge Tube	1	010-0167	100 mL Oak Ridge Tube	1	010-1114
_	6 mL Oak Ridge Tube	1	010-0382	50 mL Conical Tube	1	010-1102
	5 mL Conical Tube	1	010-1123	50 mL Oak Ridge Tube	2	010-1112
	0 mL Oak Ridge Tube	1	010-1306	30 mL Oak Ridge Tube	3	010-1115
1	0 mL BD Vacutainer™ Tube	1	010-1068	16 mL Oak Ridge Tube	7	010-1105
3	mL BD Vacutainer Tube	1	010-1128	15 mL Conical Tube	6	010-1099
1	mL BD Microtainer [™] Tube	3	010-1127	10 mL Oak Ridge Tube	7	010-1308
				10 mL BD Vacutainer Tube	7	010-1103
	_ Conical			3 mL BD Vacutainer Tube	14	010-1137
	0 mL Oak Ridge Tube	1	010-0377	1000		
	0 mL Oak Ridge Tube	1	010-1147	1000 mL		
_	6 mL Oak Ridge Tube	1	010-0376	500 mL Oak Ridge Tube	1	010-0145
	5 mL Conical Tube	1	75100378	250 mL Conical Tube	1	010-1096
	5 mL Millipore™ Filtration De		010-1340	250 mL Oak Ridge Tube	1	010-0150
	0 mL Oak Ridge Tube	1	010-1311	175 mL Nalgene Conical Bo		010-1132
1	0 mL BD Vacutainer Tube	1	010-1124	100 mL Oak Ridge Tube	3	010-1093
				50 mL Conical Tube	5	010-0180
100 n	ηL			50 mL Oak Ridge Tube	7	010-0191
_	0 mL Oak Ridge Tube	1	010-1194	30 mL Oak Ridge Tube	7	010-1095
	0 mL Oak Ridge Tube	1	010-1273	16 mL Oak Ridge Tube	15	010-1087
	6 mL Oak Ridge Tube	1	010-1272	15 mL Conical Tube	12	010-1079
1	0 mL Oak Ridge Tube	1	010-1310	10 mL Oak Ridge Tube	18	010-1307
1	0 mL BD Vacutainer Tube	1	010-1274	10 mL BD Vacutainer Tube	18	010-1415
3	mL BD Vacutainer Tube	3	010-1126	6 mL BD Vacutainer Tube	22	010-1416
1	mL BD Microtainer Tube	6	010-1125	4 mL BD Vacutainer Tube	19	010-1418
100 -				2 mL Filtration Tube	10	
	nL Ultraspeed			and 1.5 mL Conical Tube	12	010-1417
_	9 mL Tube	1	010-0189	1.8-2.7 mL BD Vacutainer T	Tube 30	010-1419
1,	3.5 mL Tube	1	010-0191	H3-LV Rotor		
250 m	h			Promega [™] Slicprep [™] 96 De	vice (/ ner run)?	018-029032
	00 mL Oak Ridge Tube	1	010-1119	Standard Microplates (28 pe		018-029031
	0 mL Conical Tube	1	75100136	2 mL Deep-well Microplates		018-029031
	0 mL Oak Ridge Tube	1	010-0138			010 020001
	0 mL Oak Ridge Tube	2	010-0130			
	6 mL Oak Ridge Tube	5	010-1072	ADAPTERS SOLD IN SETS OF 2		
	5 mL Corning [™] Conical	5	75101073	7 Actual fill volumes may vary from	nominal volume.	
	5 mL Conical Tube	5	010-1410		and the second se	
	0 mL Oak Ridge Tube	7	010-1309			
	0 mL BD Vacutainer Tube	7	010-1117			
	mL BD Vacutainer Tube	10	010-1138			
3		10	010-1100			M
				"CIENTIFIC	1.50	
					EL	

Centrifuge Rotor Maintenance

Centrifuge rotor maintenance is critical to the protection of your

samples. Leveraging more than 100 years of experience and leadership in centrifugation, our Thermo Scientific Rotor Safety Program, featuring on-site rotor inspection and safety clinics, ensures the longevity of your investment and the safety of your workplace by preventing premature rotor failure.

Thermo Scientific product representatives will evaluate the safety of your rotors and provide a comprehensive report for each rotor examined. As part of the inspection, our representatives will present information on proper rotor care and offer recommendations based upon the current rotor condition to maximize the performance of your centrifuge.

Please contact your sales representative to schedule a clinic or visit **www.thermoscientific.com/rotorsafety.**



thermoscientific.com/fiberlite

© 2010, 2011, 2012, 2013 Thermo Fisher Scientific Inc. All rights reserved. HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. Qiagen is a registered trademark of Qiagen. Promega and Slicprep are registered trademarks of Promega Corp. BD Vacutainer and BD Microtainer are registered trademarks of BD Biosciences. Millipore is a registered trademark of Millipore Corp. Beckman and Avanti are registered trademarks of Beckman, Coulter Inc. Hitachi is a registered trademark of Nissei Sangyo America. Corning is a registered trademark of Corning. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300 Austria +43 1 801 40 0 Belgium +32 53 73 42 41 China +86 21 6865 4588 or +86 10 8419 3588 France +33 2 2803 2180 Germany national toll free 0800 1 536 376 Germany international +49 6184 90 6000 India toll free 1800 22 8374 India +91 22 6716 2200 Italy +39 02 95059 554 Japan +81 45 453 9220 Netherlands +31 76 579 55 55 New Zealand +64 9 980 6700 Nordic/Baltic/CIS countries +358 9 329 10200 Russia +7 812 703 42 15 Spain/Portugal +34 93 223 09 18 Switzerland +41 44 454 12 22 UK/Ireland +44 870 609 9203 USA/Canada +1 866 984 3766

Other Asian countries +852 2885 4613 **Countries not listed** +49 6184 90 6000



Part of Thermo Fisher Scientific