



1 Crimp Neck ND8

		Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Screw Neck ND8	Screw Neck ND8	Short Thread ND9	Screw Neck ND10	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11	
		10152272 11591364	11561364	11531374 11541374	10672543 11521374	11541364 11571364	11551374	11565874 11767577 11595874 11732428	11515884	11575884 10080952 11515894 11833912 10572064 11834740 10145714 13405529 12951011 15166425 15136425 15176425 11871653 11595914 11836083 11707597 11585964 11727597 11515974 10509482 11545974	11707577 11742448 11531474	11816912 10326042 11535884 11732418 11565894 10678005 11742418	11555884	11505884	11874890
Agilent	1050								X		X				
Agilent	1050 (34 Pos. Tray)				X				X		X				
Agilent	1090										X				
Agilent	1090 (34 Pos. Tray)				X						X				
Agilent	1100								X		X				
Agilent	1200								X		X				
Agilent	1260 Infinity								X		X				
Agilent	1290 Infinity								X		X				
Agilent	G1888A														
Agilent	7673A			X					X		X				
Agilent	7683A			X			X		X		X				
Agilent	HS7694						X								
Agilent	7695A														
Agilent	79855(A)								X		X				
Agilent	5880								X		X				
Agilent	5890								X		X				
Agilent	6850 (27 Pos. Tray)								X		X				
Agilent	6850 (22 Pos. Tray)								X		X				
Agilent	6890								X		X				
Agilent	CTC HTS+HTC PAL	X			X	X			X		X				
Agilent	CTC GC PAL	X			X	X			X		X				
Agilent	CTC Combi PAL								X		X				
Agilent	Tekmar SOLAtek72														
Agilent	Archon Plugge + Trap														
Agilent	AQUAtek 70														
Agilent	7693A						X				X				
Agilent	HS7694														
Agilent	7697A														
Analytik Jena	multi N/C 3000 (TOC)														
Antec Leyden	Alexys						X		X		X				
Antec Leyden	AS 100						X		X	X	X				
Antek	736 Unisampler						X		X	X	X				
Antek	738						X		X	X	X				
Atas	GL Focus						X		X	X	X				
Beckman	501			X			X		X	X	X				
Beckman	502/502e			X			X		X	X	X				
Beckman	504	X				X			X	X	X				
Beckman	507/507e			X	X		X	X	X	X	X				
Beckman	508 (System Gold)								X						
Beckman	Marathon						X	X	X	X	X		X		
Beckman	Promis						X	X	X	X	X		X		
Beckman	Triathlon Standard Tray						X	X	X	X	X		X	X	
Beckman	Triathlon LSV Tray				X										
Beckman	Triathlon Super-LSV Tray														
Beckman	Triathlon Micro-Tray	X				X									
Bruker	LC51														
Cambridge Scientific Instruments	205 Series						X		X	X	X				
Cambridge Scientific Instruments	300 Series						X		X	X	X				
Carlo Erba	AS100	X			X	X	X	X	X	X	X		X		
Carlo Erba	AS200	X			X	X	X	X	X	X	X		X		
Carlo Erba	A200LC	X			X	X	X	X	X	X	X		X		
Carlo Erba	AS200S	X			X	X	X	X	X	X	X		X		
Carlo Erba	AS300	X			X	X	X	X	X	X	X		X		
Carlo Erba	AS800, 42 vial tray						X		X	X	X		X		
Carlo Erba	AS800, 60 vial tray		X	X	X		X		X	X	X		X	X	
Carlo Erba	HS250														
Carlo Erba	HS500														
Carlo Erba	HS800														
Carlo Erba	HS850														
Cecil Instruments	CE4800						X		X	X	X				
Cecil Instruments	AutoQuest						X	X	X	X	X	X			
CTC (LEAP)	LC PAL (216 Pos.)						X	X	X	X	X		X		
CTC (LEAP)	HTX PAL, HTC PAL, HTS PAL (200 Pos. Tray)	X					X	X	X	X	X		X		
CTC (LEAP)	HTX PAL, HTC PAL, HTS PAL (54/98 Pos. Tray)				X		X	X	X	X	X		X		
CTC (LEAP)	HTX PAL, HTC PAL, HTS PAL (32 Pos. Tray)														
CTC (LEAP)	Combi PAL (200 Pos. Tray), GC PAL (200Pos. Tray)	X													
CTC (LEAP)	Combi PAL (98 Pos. Tray), GC PAL (98 Pos. Tray)				X				X		X		X		
CTC (LEAP)	Combi PAL SPME Mode (98 Pos. Tray)								X		X		X		
CTC (LEAP)	Combi PAL (32 Pos. Tray), GC PAL (32 Pos. Tray)								X		X		X		
CTC (LEAP)	Combi PAL SPME Mode (32 Pos. Tray)														
CTC	PAL HPLC-Systems	X					X		X	X	X				
CTC	PAL Combi-xt Liquid Mode	X					X		X	X	X				
CTC	Combi-xt Headspace Option														
CTC	GC-xt Headspace Option	X					X		X	X	X				
CTC	PAL HTC-xt	X					X		X	X	X				
CTC	HTS-xt	X					X		X	X	X				
CTC	HTX-xt	X					X		X	X	X				
CTC	Combi-xt SPME Options														
CTC	A200S	X			X		X		X	X	X				
CTC	A200 LC	X			X		X	X	X	X	X		X		
CTC	HS 500														
DANI	ALS 39.80								X		X				
DANI	ALS 86.80								X		X				



3 Screw Neck ND8

		Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Screw Neck ND8	Screw Neck ND8	Short Thread ND9	Screw Neck ND10	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11
		10152272 11591364	11561364	11531374 11541374	11717567 11722408	11792408 11571364	11551374 11551374	11565874 11767577 11595874 11732428	11515884	10162512 11575884 10080952 11515894 11833912 10572064 11834740 10145714 13405529 12951011 15166425 15136425 15176425 11871653 11595914 11836083 11707597 11585964 11727597 11515974 10509482 11545974	11707577 11742448 11531474	11816912 10326042 11535884 11732418 11565894 10678005 11742418	11555884	11505884	11874890
DANI	ALS 1000									X		X			
DANI	HS38.50														
DANI	HS86.50														
DANI	Master AS									X		X			
Dimatec	Dimatoc 2000									X					
Dionex	Gina 50		X							X		X			
Dionex	AS 50			X				X		X	X				
Dionex	Summit ASI 100, Micro-Tray (192 Pos.)		X												
Dionex	Summit ASI 100, Analytical-Tray (117 Pos.)							X		X		X			
Dionex	Summit ASI 100, Semiprep.- Tray (63 Pos.)														
Dionex	Famos (LC Packings/Dionex)							X		X	X	X			
Dionex	UltraMate Analytical, cylindrical, WPS-3000 SL, 120 Pos. Rack (2ml)							X		X	X	X			
Dionex	UltraMate Analytical, conical, WPS-3000 SL, 120 (3x40) Pos. Rack (1.1ml-2ml w. Inserts)													X	
Dionex	UltraMate Micro, conical, WPS-3000 SL, 120 (3x40) Pos. Rack (250µl)				X										
Dionex	UltraMate Semipreparative, WPS-3000 SL, 66 (3x22) Pos. Rack (4ml)														
Dionex	UltraMate Nano/Cap/Micro, WPS-3000 SL, 216 (3x72) Pos. Rack (1.2ml)			X											
Dionex	ASE 200														
Dionex	ASE 40							X							
Dionex	HS-HV														
Dimatec	Dimatoc 200														
Dimatec	Dimatoc 300														
Dimatec	Dimatoc 400														
Finnigan	A200S	X				X		X		X		X			
Fisons	AS100	X		X		X		X	X	X		X			X
Fisons	AS200	X		X		X		X	X	X		X			X
Fisons	A200LC	X		X		X		X	X	X		X			X
Fisons	AS200S	X		X		X		X	X	X		X			X
Fisons	AS300	X		X	X	X		X	X	X		X			X
Fisons	AS800, 42 vial tray							X		X		X			X
Fisons	AS800, 60 vial tray		X	X	X			X		X		X			X
Fisons	HS250														
Fisons	HS500														
Fisons	HS800														
Fisons	HS850														
GE Instruments	Stevens® 900														
Gerstel	MPS2	X			X	X						X			X
Gilson	201/202							X		X					
Gilson	221/222							X		X					
Gilson	231/401							X		X					
Gilson	232/402							X		X					
Gilson	Aspec							X		X					
Gilson	Aspec Xli							X		X					
Gilson	Aspec XL4							X		X					
Gilson	221XL/222XL	X			X (only f. 221XL)	X									
Gilson	223	X				X									
Gilson	231XL/232XL/233XL	X			X (only f. 231XL)	X									
Gilson	Nano Injektor							X		X					
Gilson	235/235P/SP 235/SP 235P	X				X		X		X					
Gynkotek	Gina 50			X						X		X			
Hach Lange	IL 550 TOC-1N														
HTA	HT200H														
HTA	HT250D							X		X	X	X			X
HTA	HT2801							X		X	X	X			X
HTA	HT300A							X		X	X	X			X
HTA	HT310A							X		X	X	X			X
HTA	HT300L							X		X	X	X			X
ICI	LC1600	X				X									
IMT GmbH	VSP4000														
IMT GmbH	PTA3000														
Jasco	AS 2055/AS 2055 (i)				X			X		X	X	X			
Jasco	AS 2057/AS 2057 (i)				X			X		X	X	X			
Jasco	AS 2059				X			X		X	X	X			
Knauer	K-3800 (Basic Marathon)							X		X	X	X			X
Knauer	Smartline K-3950							X		X	X	X			
Knauer	PLATInblue AS-1							X		X	X	X			
LDC	713-60		X		X										
LDC	Marathon							X	X	X		X			X
LDC	Promis							X	X	X		X			X
LEAP	pls. see CTC														
O.I. Analytical	1020A														
O.I. Analytical	1088														
O.I. Analytical	1096+														
O.I. Analytical	4551A														
O.I. Analytical	1552														
PerkinElmer	Series 200, 25 vial tray											X			X
PerkinElmer	Series 200, 85 vial tray											X			X
PerkinElmer	Series 200, 81/100 vial tray											X			X



7 Crimp Neck ND11

		Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Screw Neck ND8	Screw Neck ND8	Short Thread ND9	Screw Neck ND10	Crimp Neck ND11 7	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11
		10152272 11591364	11561364	11531374 11541374	10672543 11521374	11541364 11571364	11551374	11565874 10560053 11595874 11525884	11515884	10162512 11575884 10080952 11515894 11575894 11585894 11595894 10145714 11535914 12951011 15166425 15136425 15176425 11575914 11595914 11515924 11565964 11585964 11505974 11515974 10509482 11545974	11511474 10521593 11531474	10081022 10326042 11535884 11545884 11565894 10678005 11505894	11555884 11555894	11505884	11555894
Spark	Midas, Standard 84 Pos. Tray							X	X	X		X		X	X
Spark	Midas, Large Capacity 96 Pos. Tray							X	X	X		X		X	X
Spark	Midas, Large Volume 24 Pos. Tray														
Spark	Alias							X	X	X		X		X	
Spark	Promis							X	X	X		X		X	
Spark	SPH 125							X	X	X		X		X	
Spark	Triathlon, Standard 96 Tray							X	X	X		X		X	X
Spark	Triathlon, LSV 72 Pos. Tray														
Spark	Triathlon, Super-LSV 32 Pos. Tray														
Spark	Triathlon, Micro 160 Pos. Tray	X				X									
Spark	Endurance 48 Pos. Tray							X	X	X		X		X	
Spark	Reliance 48 Pos. Tray							X	X	X		X		X	
Spark	Integrity									X		X			
Spark	Prospekt 2											X			
Spark	Reliance/Symbiosis Pharma											X			
Spark	Dried Blood Spot (DBS)							X		X		X			
Spark	Integrity							X		X		X			
Spark	Integrity 108 Pos.(2ml) 2 x Plates , IntegrityPlus 2 x 108 Pos (2ml) 4 x Plates							X		X		X			
Spark	Optimas							X		X		X			
Spark	Optimas 96 Pos. (2mL) 24 Pos. (10mL)							X		X		X			
Spark	Alias							X	X	X		X		X	
Spectra-Physics	8875							X		X		X			
Spectra-Physics	8880							X		X		X			
Spectra-Physics	SpectraSYSTEM AS1000	X			X	X		X		X		X		X	
Spectra-Physics	SpectraSYSTEM AS 3000	X		X	X	X		X		X		X		X	X
Spectra-Physics	SpectraSYSTEM AS 3500	X			X	X		X		X		X		X	
Sykam	S 5200									X					
Teledyne Tekmar	7000/7000HT/7050														
Teledyne Tekmar	AQUATEK 70/SOLATEK 72™														
Teledyne Tekmar	STS 8000 TOC														
Teledyne Tekmar	HT3														
Thermo Scientific	AS1000 (Trace GC)	X			X	X		X		X		X		X	
Thermo Scientific	AS200	X			X	X		X	X	X		X		X	
Thermo Scientific	AS300	X			X	X		X	X	X		X		X	
Thermo Scientific	AS2000 30 vial tray														
Thermo Scientific	AS2000 90 vial tray (Trace GC)			X				X	X	X		X		X	X
Thermo Scientific	A13000 (II)/AS3000 (II) AS3500 (Trace GC + Focus GC)	X			X	X				X		X		X	
Thermo Scientific	A200LC	X				X		X	X	X		X		X	
Thermo Scientific	SpectraSYSTEM AS 1000	X			X	X		X		X		X		X	
Thermo Scientific	SpectraSYSTEM AS 3000	X		X	X	X		X		X		X		X	X
Thermo Scientific	SpectraSYSTEM AS 3500	X			X	X		X		X		X		X	
Thermo Scientific	A200S	X			X	X		X	X	X		X		X	
Thermo Scientific	AS100	X			X	X		X	X	X		X		X	
Thermo Scientific	AS800, 42 vial tray							X		X		X		X	X
Thermo Scientific	AS800, 60 vial tray							X		X		X		X	
Thermo Scientific	Dionex AS-AP		X	X	X			X		X		X			
Thermo Scientific	Dionex UltiMate WPS-3000	X		X				X	X	X		X			
Thermo Scientific	Dionex AS 40														
Thermo Scientific	HS250														
Thermo Scientific	HS500														
Thermo Scientific	HS800														
Thermo Scientific	HS850														
Thermo Scientific	HS2000														
Thermo Scientific	TriPlus (=GC PAL) (AS+ Duo)	X		X	X	X		X		X		X		X	X
Thermo Scientific	TriPlus HS														
Thermo Scientific	TriPlus SPME														
Thermo Scientific	TriPlus RSH				X	X		X		X		X			
Thermo Scientific	TriPlus 300														
Thermo Scientific	HiPerTOC														
Thermo Scientific	Surveyor (Surveyor Plus)	X			X	X		X		X		X		X	
Thermo Scientific	Accela High Speed LC Autosampler (200 Pos.)	X			X	X		X		X		X			
Thermo Scientific	Accela Open Autosampler	X				X	X		X	X		X			
Thermo Scientific	Trace 1300 Series							X		X	X	X			
Unicam	4247							X		X		X			
Unicam	4710							X		X		X			
Unicam	4700 (GC)	X				X									
Unicam	4700 (LC)	X				X		X	X	X					
Unicam	LC-XP							X	X	X		X		X	
Unicam	S4/S8	X				X									
Varian	ProStar 400, Standard 96 Pos. Tray							X	X	X	X	X		X	
Varian	ProStar 400, King Size 48 Pos. Tray														
Varian	ProStar 410, Standard 84 Pos. Tray							X	X	X	X	X		X	X
Varian	ProStar 410, Large Capacity 96 Pos. Tray							X	X	X	X	X		X	X
Varian	ProStar 410, Large Volume 24 Pos. Tray														
Varian	ProStar 420, Standard 96 Pos. Tray							X	X	X	X	X		X	X



9 Short Thread ND9

		Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Crimp Neck ND8	Screw Neck ND8	Screw Neck ND8	Short Thread ND9	Screw Neck ND10	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11	Crimp Neck ND11	
		10152272 11591364	11561364	11531374 11541374	10672543 11521374	11541364 11571364	11551374	11565874 10560053 11595874 11525884	11515884	10162512 11575884 10080952 11515894 11575894 11585894 11595894 10145714 11535914 12951011 15166425 15136425 15176425 11575914 11595914 11515924 11565964 11585964 11505974 11515974 10509482 11545974	11511474 10521593 11531474	10081022 10326042 11535884 11545884 11565894 10678005 11505894	11555884	11505884	11555894	
Varian	ProStar 420, LSV 72 Pos. Tray				X											
Varian	ProStar 420, Super-LSV 32 Pos. Tray															
Varian	ProStar 420, Micro 160 Pos. Tray	X				X										
Varian	ProStar 430, 48 Pos. Tray							X	X	X		X			X	
Varian	8035							X	X	X						
Varian	8000							X		X		X				
Varian	8400 (100 Pos.)							X		X		X				
Varian	8410-Autoinjector (10 x 2ml; 6 x 5ml; 5 x 10ml)							X		X		X				
Varian	8100							X		X		X				
Varian	8200							X		X	X	X				
Varian	CP-910, 911, 912							X		X	X	X				
Varian	CP-940, 941											X				
Varian	LC 9100/LC 9095/LC 9090									X		X				
Varian	Archon															
Varian	COMBI PAL (200 Pos. Tray) GC PAL (200 pos. Tray)	X				X										
Varian	COMBI PAL (98 Pos. Tray) GC PAL (98 Pos. Tray)				X					X		X				
Varian	COMBI PAL SPME mode (98 Pos. Tray)									X		X				
Varian	COMBI PAL (32 Pos. Tray) GC PAL (32 Pos. Tray)															
Varian	COMBI PAL SPME mode (32 Pos. Tray)															
Varian	Genesis															
Varian	Marathon Basic, Standard 96 Pos. Tray							X	X	X		X			X	
Varian	Marathon Basic, Präp, King Size 48 Pos. Tray															
Varian	Vista							X		X						
Varian	CP-9020/CP-9025															
Varian	CP-9060															
Varian	CP-9010							X		X		X				
Varian	920-LC/940-LC							X								
Viscotek	Vortex™															
Viscotek	GPC Autosampler							X		X	X					
WWR	Hitachi Chromaster							X		X		X				
WWR(Merck®)/Hitachi	L2200 (LaChrom Elite)/L2200-U (LaChrom Ultra) (200 Pos. Tray)							X		X						
WWR(Merck®)/Hitachi	L2200 (LaChrom Elite) (128 Pos. Tray)															
WWR(Merck®)/Hitachi	L7200 (LaChrom) (80 Pos. Tray)/L7250(LaChrom) (120 Pos. Tray)							X		X						
WWR(Merck®)/Hitachi	L7250 (LaChrom) (Rack Holder for combination Racks)				X			X		X						
WWR(Merck®)/Hitachi	655-A40 (108 Pos. Tray)							X		X						
WWR(Merck®)/Hitachi	L-9100							X		X						
WWR(Merck®)/Hitachi	AS 2000 (50 Pos. Tray)							X		X		X				
WWR(Merck®)/Hitachi	AS 4000 (150 Pos. Tray)							X		X		X				
WWR(Merck®)/Hitachi	AS 4000 (198 Pos. Tray)				X											
WWR(Merck®)/Hitachi	5210 (Chromaster) 195 Pos (1 mL), 120 Pos 1.5 mL (Standard), 72 Pos. (4 mL), 2 x MTP (96,384)							X		X		X				
WWR(Merck®)/Hitachi	AS 6000				X			X		X						
Waters®	ACQUITY™ UPLC Systeme									X						
Waters®	Wisp 48 position															
Waters®	Wisp 96 position															
Waters®	717, 96 Position Carousel															
Waters®	717, 48 Position Carousel															
Waters®	Alliance®									X	X	X				
Waters®	Alliance® GPC 2000															
Waters®	Alliance® HT Syst.									X	X	X				
Waters®	Alliance® 2790/2795									X	X	X				
Waters®	Acquity Sample Organizer									X						
Waters®	Acquity/CapLC/Waters/Nano Acquity									X						
Waters®	Alliance® 2690/2695									X (not for 11515894)	X	X				



10 Snap Ring ND11

		Snap Ring ND11	Screw Neck ND13	Shell Vials	Shell Vials	Shell Vials	Shell Vials	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Screw Neck ND24
		11525894	10571013	10306062	10145424	11555914	11516074	10663303	10681033	12990951	10680843	11506114	10192652	11520545	10465982
		11535894	11556044		10506075	11565914	10455982		10195012	10080822	12971231	12941221			10000782
		11545894	11576044		11561374					10152512		11526114			10758874
		11545914	11586044		10224852					12981241		10070952			11510585
		12672465										10510323			10458082
		15219468										12910991			11530585
		11585914													11540585
		11505924													11550585
		11525924													
		11545924													
		11575964													
		11595964													
		10488082													
		10705203													
Varian	ProStar 420, LSV 72 Pos. Tray		X												
Varian	ProStar 420, Super-LSV 32 Pos. Tray							X	X						
Varian	ProStar 420, Micro 160 Pos. Tray														
Varian	ProStar 430, 48 Pos. Tray	X													
Varian	8035														
Varian	8000														
Varian	8400 (100 Pos.)												X		
Varian	8410-Autoinjector (10 x 2ml, 6 x 5ml, 5 x 10ml)												X	X	
Varian	8100														
Varian	8200														
Varian	CP-910, 911, 912														
Varian	CP-940, 941														
Varian	LC 9100/LC 9095/LC 9090	X													
Varian	Archon														X
Varian	COMBI PAL (200 Pos. Tray)								X						
Varian	GC PAL (200 pos. Tray)									X					
Varian	COMBI PAL (98 Pos. Tray)									X					
Varian	GC PAL (98 Pos. Tray)										X				
Varian	COMBI PAL SPME mode (98 Pos. Tray)									X					
Varian	COMBI PAL (32 Pos. Tray)									X			X		
Varian	GC PAL (32 Pos. Tray)											X			
Varian	COMBI PAL SPME mode (32 Pos. Tray)									X		X			
Varian	Genesis										X				
Varian	Marathon Basic, Standard 96 Pos. Tray	X													
Varian	Marathon Basic, Präp, King Size 48 Pos. Tray							X							
Varian	Vista														
Varian	CP-9020/CP-9025								X						
Varian	CP-9060								X						
Varian	CP-9010														
Varian	920-LC/940-LC	X													
Viscotek	Vortex™														X
Viscotek	GPC Autosampler														
WWR	Hitachi Chromaster		X												
WWR(Merck®)/Hitachi	L2200 (LaChrom Elite)/L2200-U (LaChrom Ultra) (200 Pos. Tray)	X													
WWR(Merck®)/Hitachi	L2200 (LaChrom Elite) (128 Pos. Tray)		X												
WWR(Merck®)/Hitachi	L7200 (LaChrom) (80 Pos. Tray)/L7250(LaChrom) (120 Pos. Tray)	X													
WWR(Merck®)/Hitachi	L7250 (LaChrom) (Rack Holder for combination Racks)	X	X												X
WWR(Merck®)/Hitachi	655-A40 (108 Pos. Tray)	X													
WWR(Merck®)/Hitachi	L-9100	X													
WWR(Merck®)/Hitachi	AS 2000 (50 Pos. Tray)	X													
WWR(Merck®)/Hitachi	AS 4000 (150 Pos. Tray)	X													
WWR(Merck®)/Hitachi	AS 4000 (198 Pos. Tray)														
WWR(Merck®)/Hitachi	5210 (Chromaster) 195 Pos (1 mL), 120 Pos 1.5 mL (Standard), 72 Pos. (4 mL), 2 x MTP (96,384)	X	X		X		X								
WWR(Merck®)/Hitachi	AS 6000	X													
Waters®	ACQUITY™ UPLC Systeme		X (for 24 Position Plate)												
Waters®	Wisp 48 position		X				X								
Waters®	Wisp 96 position				X										
Waters®	717, 96 Position Carousel				X										
Waters®	717, 48 Position Carousel		X				X								
Waters®	Alliance®	X													
Waters®	Alliance® GPC 2000		X						X						
Waters®	Alliance® HT Syst.	X													
Waters®	Alliance® 2790/2795	X													
Waters®	Acquity Sample Organizer	X													
Waters®	Acquity/CapLC/Waters/Nano Acquity	X													
Waters®	Alliance® 2690/2695	X													

According to our information our compatibility chart was correct when it was printed. Please note that later changes in manufacturer's specifications may lead to different compatibilities.

© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500



A Thermo Fisher Scientific Brand